



2021 Product List

Diagnostic Automation / Cortez Diagnostic Inc.

CLIA



DA /CD

CLIA

Product Name	Catalog #	# of Tests	CE	IVD	Qualitative / Quantitative	Assay Time	Sensitivity
Allergy Assays							
AccuDiag™ IgE	9050-16	96	✓	IVD	Quantitative	~65 min	5 IU/mL
Cancer (Tumor Marker) Assays							
AccuDiag™ AFP	9040-16	96	✓	IVD	Quantitative	~65 min	2.0 ng/mL
AccuDiag™ CEA	9041-16	96	✓	IVD	Quantitative	~65 min	1.0 ng/mL
AccuDiag™ PSA	9042-16	96		In-Vitro	Quantitative	~125 min	0.5 ng/mL
AccuDiag™ Free Beta hCG	9043-16	96	✓	IVD	Quantitative	~65 min	0.1 mIU/mL
AccuDiag™ Free PSA	9054-16	96		In-Vitro	Quantitative	~125 min	0.05 ng/mL
AccuDiag™ Beta 2- Microglobulin	9076-16	96	✓	IVD	Quantitative	~65 min	2.0 µg/mL
AccuDiag™ NSE	9056-16	96	✓	IVD	Quantitative	~65 min	2.0 ng/mL
AccuDiag™ CA-12-5	9057-16	96	✓	IVD	Quantitative	~125 min	5 U/mL
AccuDiag™ CA-19-9	9059-16	96	✓	IVD	Quantitative	~125 min	5 U/mL
AccuDiag™ CA-15-3	9075-16	96	✓	IVD	Quantitative	~125 min	5 U/mL
AccuDiag™ Ferritin	9019-16	96	✓	IVD	Quantitative	~75 min	5.0 ng/mL
Cardio-Vascular Monitoring Assays							
AccuDiag™ CRP	9053-16	96	✓	IVD	Quantitative	~125 min	5 U/mL
Diabetes							
AccuDiag™ Insulin	9021-16	96	✓	IVD	Quantitative	~65 min	2 µIU/mL
AccuDiag™ C-peptide	9022-16	96	✓	IVD	Quantitative	~65 min	0.05 ng/mL
Fertility Assays							

CLIA

Product Name	Catalog #	# of Tests	CE	IVD	Qualitative / Quantitative	Assay Time	Sensitivity
AccuDiag™ LH	9012-16	96	✓	IVD	Quantitative	~65 min	2 mIU/mL
AccuDiag™ FSH	9013-16	96	✓	IVD	Quantitative	~65 min	2.5mIU/mL
AccuDiag™ Prolactin (PRL)	9014-16	96	✓	IVD	Quantitative	~65 min	2 ng/mL
AccuDiag™ hCG	9015-16	96	✓	IVD	Quantitative	~65 min	2 mIU/mL

Infectious Disease Assays

AccuDiag™ H. pylori IgA Min. Order: 5 kits	9060-11	96	✓	IVD	Qualitative	~65 min	N/A
AccuDiag™ H. pylori IgG Min. Order: 5 kits	9061-11	96	✓	IVD	Qualitative	~65 min	99%
AccuDiag™ H. pylori IgM Min. Order: 5 kits	9062-11	96	✓	IVD	Qualitative	~65 min	100%
AccuDiag™ Toxoplasma IgG Min. Order: 5 kits	9071-11	96		In-Vitro	Qualitative	~65 min	94%
AccuDiag™ Toxoplasma IgM Min. Order: 5 kits	9072-11	96		In-Vitro	Qualitative	~65 min	100%
AccuDiag™ CMV IgG Min. Order: 5 kits	9063-11	96		In-Vitro	Quantitative	~65 min	97%
AccuDiag™ CMV IgM Min. Order: 5 kits	9064-11	96		In-Vitro	Qualitative	~65 min	100%
AccuDiag™ HSV1 IgG Min. Order: 5 kits	9065-11	96	✓	IVD	Qualitative	~65 min	100%
AccuDiag™ HSV1 IgM Min. Order: 5 kits	9066-11	96	✓	IVD	Qualitative	~65 min	100%
AccuDiag™ HSV2 IgG Min. Order: 5 kits	9067-11	96	✓	IVD	Qualitative	~65 min	100%
AccuDiag™ HSV2 IgM Min. Order: 5 kits	9068-11	96	✓	IVD	Qualitative	~65 min	70%
AccuDiag™ Rubella IgG Min. Order: 5 kits	9069-11	96		In-Vitro	Qualitative	~65 min	100%
AccuDiag™ Rubella IgM Min. Order: 5 kits	9070-11	96		In-Vitro	Qualitative	~65 min	100%

Steroids Assays

AccuDiag™ Progesterone	9025-16	96	✓	IVD	Quantitative	~95 min	0.2 ng/mL
AccuDiag™ Estradiol	9028-16	96	✓	IVD	Quantitative	~95 min	5 pg/mL
AccuDiag™ Testosterone	9052-16	96	✓	IVD	Quantitative	~95 min	0.05 ng/mL

CLIA

Product Name	Catalog #	# of Tests	CE	IVD	Qualitative / Quantitative	Assay Time	Sensitivity
AccuDiag™ Thyroid Assays							
AccuDiag™ T3	9001-16	96	✓	IVD	Quantitative	~65 min	0.03 ng/mL
AccuDiag™ T4	9003-16	96	✓	IVD	Quantitative	~50 min	0.5 µg/mL
AccuDiag™ TSH	9005-16	96	✓	IVD	Quantitative	~65 min	0.03 ng/mL